# Silosept®

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\*Selon les recommandations  
RKI- Bundesgesundheitsblatt 01/2004

## Immersion disinfectant for dental impressions and restorations.

**Range of indications:** Immersing silicone, vinylsiloxanether, polyether and alginate dental impressions as well as dental restorations to disinfect.

**Properties:** The mixed solution is microbiologically highly effective and cleans well in only 10 minutes immersion time – it is compatible with impression materials, which remain dimensionally stable.

**Mode of application:** Sprinkle 2 scoop (20 g) of Silosept® into 2 litres of lukewarm water (creates a 2 % solution). Rinse the impressions and restorations in running water and immerse them in the solution. Allow 10 minutes for it to take effect, remove them and rinse again in running water.

### Change disinfection solution daily.

**Microbiological effectiveness:** A 2 % Silosept® solution deactivates the relevant germ spectrum – bacteria (incl. Tb) and fungi, enveloped viruses\* (incl. HIV, HBV, HCV) – in 10 minutes, antiviral agent at 2 % - 30 minutes and is easily biodegradable.

**Special instructions:** If the patient has been diagnosed with TB, extend the reaction time to 30 minutes and use impression materials that are appropriate for this application. Test material compatibility prior to use. Direct contact with the skin and fabrics may result in bleaching effects. The possibility of sensitive individuals experiencing allergic reactions to the product cannot be ruled out. Wear protective gloves when placing/removing the impressions and restorations.

**Composition:** 100 g Silosept® contain: 45 g pentapotassium-bis(peroxymonosulfate)bis(sulphate)

**Labelling according to Regulation (EC) No. 648/2004:** 5 – 15 % anionic surfactants, < 5 % non-ionic surfactants, < 5 % soap and phosphonates, perfume.

**Please note: R36/38** - Irritating to eyes and skin. **R52** Harmful to aquatic organisms. **S26** - In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. **S28** After contact with skin, wash immediately with plenty of water.

Do not store above room temperature and keep closed. Do not use after expiration date.

\*In accordance with the Rôbert Koch- Institute guideline - German Federal Health Gazette 01/2004

CE 0123

800 g

REF 18701

Made in  
Germany